

# Application Data Sheet

## No. 180

### System Gas Chromatograph

#### Octans and Lower Boiling Hydrocarbons in Olefin-Free Gasolines

Nexis GC-2030LBH3

This system is for determining the composition of C8 or lower boiling paraffins and naphthenes in olefin-free gasolines. Detailed identification is enabled by comparing chromatograms of two different analytical conditions.

#### Analyzer Information

##### System Configuration:

One SPL injector / one capillary column / one FID

##### Methods met:

UOP-690

#### System Features

- Single FID channel
- Good repeatability
- Detailed hydrocarbon analysis with simple configuration
- Dedicated software for supporting identification

#### Typical Chromatograms




Fig. 1 Chromatogram of FID - 1 of 4




Fig. 2 Chromatogram of FID - 2 of 4




Fig. 3 Chromatogram of FID - 3 of 4




Fig. 4 Chromatogram of FID - 4 of 4

ID#	Name	ID#	Name
1	Isobutane	30	2,5-Dimethylhexane
2	n-Butane	31	2,2,3-Trimethylpentane + 2,4-Dimethylhexane
3	Isopentane	32	1-trans-2-cis-4-Trimethylcyclopentane
4	n-Pentane	33	3,3-Dimethylhexane
5	2,2-Dimethylbutane	34	1-trans-2-cis-3-Trimethylcyclopentane
6	Cyclopentane	35	2,3,4-Trimethylpentane
7	2,3-Dimethylbutane	36	2,3,3-Trimethylpentane + Toluene
8	2-Methylpentane	37	2,3-Dimethylhexane + 1,1,2-Trimethylcyclopentane
9	3-Methylpentane	38	2-Methyl-3-ethylpentane
10	n-Hexane	39	2-Methylheptane
11	2,2-Dimethylpentane	40	4-Methylheptane
12	Methylcyclopentane	41	3,4-Dimethylhexane + 3-Methyl-3-ethylpentane
13	2,4-Dimethylpentane	42	3-Methylheptane
14	2,2,3-Trimethylbutane	43	3-Ethylhexane + 1-cis-3-Dimethylcyclohexane + 1-cis-2-trans-3-Trimethylcyclopentane
15	Benzene	44	1-trans-4-Dimethylcyclohexane
16	3,3-Dimethylpentane	45	1,1-Dimethylcyclohexane
17	Cyclohexane	46	1-Methyl-trans-3-ethylcyclopentane
18	2-Methylhexane	47	1-Methyl-cis-3-ethylcyclopentane
19	2,3-Dimethylpentane	48	1-Methyl-trans-2-ethylcyclopentane
20	1,1-Dimethylcyclopentane	49	1-Methyl-1-ethylcyclopentane
21	3-Methylhexane	50	1-trans-2-Dimethylcyclohexane
22	1-trans-3-Dimethylcyclopentane	51	n-Octane + 1-cis-4-Dimethylcyclohexane + 1-trans-3-Dimethylcyclohexane
23	3-Ethylpentane	52	Isopropylcyclopentane
24	1-trans-2-Dimethylcyclopentane	53	1-Methyl-cis-2-ethylcyclopentane
25	2,2,4-Trimethylpentane	54	1-cis-2-Dimethylcyclohexane
26	n-Heptane	55	n-Propylcyclopentane + Ethylcyclohexane
27	Methylcyclohexane + 1-cis-2-Dimethylcyclopentane	56	Ethylbenzene
28	1,1,3-Trimethylcyclopentane + 2,2-Dimethylhexane	57	m-Xylene + p-Xylene
29	Ethylcyclopentane	58	o-Xylene

Fig. 5 Compound List of the Chromatogram

First Edition: November,2017



For Research Use Only. Not for use in diagnostic procedures.  
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.