

Hydrogen/Carbon Dioxide Monitoring System

Introduction


Advancements in H₂ production, carbon capture and CO₂ recycling techniques require robust measurement and analysis.

Electrocatalysis, photocatalysis and other methods produce key products including H₂, O₂, N₂, CH₄, CO, CO₂, and valuable hydrocarbons.

Shimadzu provides a Gas Chromatography system that facilitates easy qualitative and quantitative analysis of these major components.



Nexis GC-2030 with 2 detectors & top-mounted valves



*The chromatograms shown are based on example data. *Other configurations can affect results.*

System Overview


Shimadzu's flagship GC-2030 configured for H₂/CO₂ Reactor Monitoring offers:

- Bespoke build
- Sample run-time less than 10 minutes
- Fully developed analytical method
- Flexible and future-proof options with upgrade pathways available
- SPL injection port, FID detector for valuable light hydrocarbons C₂-C₅
- TCD detector and 10-port valve for permanent gas analysis
- Installation by factory-trained Shimadzu Engineer
- On-site bespoke training by Shimadzu Application Specialist
- Ongoing technical support online for the instrument's lifetime
- Windows 10/11 PC with intuitive LabSolutions Software package

System Schematic

Using a 10-port valve (V1) with backflush capability, using Column A with TCD detection for the analysis of permanent gases including H₂, O₂, N₂, CH₄, CO and CO₂. Expected LODs = 100 ppm (lower achievable, dependent on carrier gas used).

A 6-port valve (V2) injects onto Column B with FID detection for the analysis of light hydrocarbons (chain length C₂ to C₅). Expected LODs = 20 ppm.



Flexible Options and Future-Proofing

- A** **Column options depending on analytical requirements**
Column A choice: Molsieve 5A for baseline separation of all permanent gases, or Carboxen 1010 to facilitate separation of H₂, O₂, CO, CH₄ and CO₂ on the TCD.
- B** **High-sensitivity analysis with BID**
Use BID (Barrier Ionisation Discharge) detector to improve LODs of H₂, O₂, N₂, CH₄, CO and CO₂ by replacing the TCD.
- C** **Liquid Sample analysis (valuable alcohols and hydrocarbons)**
Column B switched out for hydrocarbon capillary column for analysis of liquid samples.
- D** **CO₂ detection using methanizing technology**
Replace the FID nozzle with Shimadzu's Jetanizer™ to facilitate CO₂ detection on FID.



Shimadzu UK Limited

Email: info@shimadzu.co.uk Phone: 01908 552209 Web: www.shimadzu.co.uk