

Prices valid until 1<sup>st</sup> January 2010.

## Syringes

Part No.	Shimadzu No.	Description	Price (£)	Qty
0156000	980-02895	0.5µL (plunger-in needle) syringe with 0.63mm removable needle	60.00	
0156005*	980-02896	5µL syringe with 0.63mm fixed needle	44.00	
0156011*	980-02878	10µL syringe with 0.63mm fixed needle	40.00	
0156013	980-02898	10µL syringe with 0.63/0.47mm tapered fixed needle	40.00	
0156015	980-02899	10µL Gas Tight syringe with 0.63mm removable needle	44.00	
0156300	980-02900	Replacement plunger and needle kit for P/N 015600 (pkt 1)	48.00	
0156315	980-02901	Replacement needle for P/N 0156015 (pkt 2)	16.00	
0156316	980-02902	Replacement plunger for P/N 0156015 (pkt 2)	26.00	

## GC, GC-MS

Part No.	Shimadzu No.	Description	Price (£)	Qty
041890SH*	980-02879	Enduro Blue Septum (pkt 50)	72.00	
072663SH*	980-02935	15% Graphite/85% Vespel ferrule for 0.1-0.25mm ID columns (pkt 10)	29.00	
072654SH*	980-02932	15% Graphite/85% Vespel ferrule for 0.32mm ID columns (pkt 10)	29.00	
072655SH*	980-02933	15% Graphite/85% Vespel ferrule for 0.53mm ID columns (pkt 10)	29.00	
0726080SH*	980-02928	100% Graphite ferrule for 0.1-0.32mm ID columns (pkt 10)	70.00	
0726082SH*	980-02929	100% Graphite ferrule for 0.45-0.53mm ID columns (pkt 10)	70.00	
0726001SH*	980-04039	100% Graphite ferrule for 5mm ID packed columns (pkt 10)	34.00	
073200SH*	980-02858	Siltite Metal ferrule for 0.1-0.25mm ID columns (2x nuts, 10x ferrules)	53.00	
073201SH	980-04040	Siltite Metal ferrule for 0.32mm ID columns (2x nuts, 10x ferrules)	53.00	
073202SH	980-04041	Siltite Metal ferrule for 0.53mm ID columns (2x nuts, 10x ferrules)	53.00	
073210SH	980-04042	Siltite Metal ferrule for 1/32" ID columns (2x nuts, 10x ferrules)	53.00	
073220SH	980-02859	Replacement Siltite Metal ferrule for 0.1-0.25mm ID columns (pkt 10)	34.00	
073221SH	980-02938	Replacement Siltite Metal ferrule for 0.32mm ID columns (pkt 10)	34.00	
073222SH	980-02939	Replacement Siltite Metal ferrule for 0.53mm ID columns (pkt 10)	34.00	
073219SH	980-04043	Replacement Siltite Metal ferrule for 1/32" ID columns (pkt 10)	34.00	
073224SH	980-02860	Replacement Siltite Metal nut (pkt 5)	48.00	

**GC, GC-MS cont.**

Part No.	Shimadzu No.	Description	Instrument	Price (£)	Qty
092064SH*	980-04045	Split liner, straight-through - 3.4mm ID (pkt 5)	2014, 2010, 17A	57.00	
092061SH	980-04046	Split / splitless liner with gooseneck & quartz wool - 3.4mm ID (pkt 5)	2014, 17A	94.00	
092071SH	980-02944	Split / splitless liner with single taper - 3.4mm ID (pkt 5)	2014, 2010, 17A	127.00	
092077SH	980-02945	Split / splitless liner with middle gooseneck - 3.4mm ID (pkt 5)	2010	60.00	
092062SH	980-02942	Straight-through Focusliner - 3.4mm ID (pkt 5)	2014, 17A	97.00	
092068SH	980-02880	Tapered Focusliner - 3.4mm ID (pkt 5)	2014, 17A	133.00	
092085SH	980-02946	Split / splitless liner with middle gooseneck - 3.4mm ID (pkt 5)	2014, 17A	94.00	
0920861SH*	980-02947	Splitless liner, straight-through - 2.6mm ID (pkt 5)	2014, 2010, 17A	50.00	
092087SH	980-04047	Direct injection liner for 0.53mm ID columns - 2.6mm ID (pkt 5)	2014, 17A	125.00	
092089SH	980-04048	SPME liner - 0.75mm ID (pkt 5)	2014, 17A	68.00	
0726001SH	980-04039	Graphite Sealing Ring (pkt 10)	SPL-14	34.00	
0726007SH	980-02927	Graphite Sealing Ring (pkt 10)	SPL-17	35.00	
0726533SH*	980-02931	Viton O-Ring (pkt 10)	SPL-2014 & 2010	19.00	
054101SH*	980-02830	SHIM-5MS Capillary Column 30m x 0.25mm ID x 0.25µm		347.00	

\* included in standard kit

**Other**

Part No.	Shimadzu No.	Description	Price (£)	Quantity

**For GC Capillary Columns pricing please contact Shimadzu UK Limited**

Customer Ref No..... Department.....

Company / Organisation..... Contact Name.....

..... Order.....

Name..... Signature.....